



## Human AB Serum

91-02-0023 100ml

biological source : from human male AB plasma

sterility sterile-filtered

composition hemoglobin,  $\leq 20\text{mg/dL}$

origin USA origin

impurities  $\leq 10\text{ EU/mL}$  endotoxin

shipped in dry ice

storage temp.  $-20^{\circ}\text{C}$

Analysis Note Endotoxin tested

### Warning

Each donor is tested for and found non-reactive for Hepatitis B & C and non-reactive for Human Immunodeficiency Virus(HIV) antibody by ELISA. Nevertheless, products of human origin should be considered potentially infectious and handled accordingly.

### Application

Human AB serum is used in tissue engineering, transplantation and cell therapy applications for the expansion of mesenchymal stem cells (MSC) from adipose tissue or mesenchymal stromal cells from human bone marrow; for standardized limbal epithelial stem cell graft generation and transplantation; for ex vivo expansion of NK cells from peripheral blood in Hematopoietic Stem Cell Expansion Medium; and for upgrading pretransplant human islet culture technology.

苏州欣协生物科技有限公司

电话: 0512-63037851

网址: [www.xinbio.cn](http://www.xinbio.cn)

邮箱: [info@szxxbio.com](mailto:info@szxxbio.com)